



Uricult[®] Vet


Quick Procedure Guide

1. Collect specimen by cystocentesis or catheterization.
2. Unscrew cap of Uricult Vet vial and place 2 drops of urine from the syringe on each side of the media paddle. Rotate the paddle to be sure media surface is completely covered.
3. Re-Insert the Uricult Vet dip slide into its plastic vial and loosely close the cap.
4. Place the vial into an incubator set at 36 +/- 2° C for 16 - 24 hours.
5. Remove from incubator and compare both sides of media to colony density chart and check for color change.
6. Use the Uricult Vet Slide Guide for preliminary (presumptive) identification and gram (+) or gram (-) results. (Slide - Guide available from LifeSign, LLC)

Read the package insert completely before performing test

Visit www.uricultusa.com for training video and additional information

INTERPRETATION OF RESULTS Colony Density Chart



Expected values

Method of sampling	Significant colony count CFU/ml	
	Dog	Cat
Bladder aspiration	≥ 1,000	≥ 1,000
Catheterization	≥ 10,000	≥ 1,000

For more information about Uricult[®] Vet or to place an order, contact:



85 Orchard Road, Skillman, New Jersey 08558
800.526.2125 732.246.3366 Fax: 732.246.0570
www.lifesignmed.com